

Quick Manual for Asalab and equipment for ECoG

This manual is not a complete description of ANT equipment for recording and Asalab software. It was made in order to help for a quick use with basic knowledge of the software and equipment. For more details there is the manual of the software in the help tab in Asalab and for the equipment there is all the ANT documentation in the laboratory.

1) Equipment

a- Amplifier and cables


In Figure 1, here is the amplifier with all the cables to connect it to a computer:



This is the FUSBI, it's an interface between glass fiber and USB port:




This orange cable is the glass fiber. The must not be bent or folded!!!



The blue box is the amplifier power supply box.




This blue box is connected to the amplifier by the grey cable.



b- Filter

It's an high-pass filter in front of the normal analogue input. We will use this filter because the electrode on ECoG grids are made from Platina, whereas other metals on the body may be used for example Ag/AgCl. And it may cause drifts on the signal.

It is an RC filter (resistance capacity) and its cut-off frequency is 1.6 Hz. For more details, see Annex 1.



In case of recording only 32 channels, Input 2 and Output 2 might be unused. In this case it is recommended to short circuit all pins on the input connector of the filter bank. The white box function is to short circuit all these empty connectors. (The output 2 connector must not be connected to the amplifier in this case).

In the next pictures, we can see the connector where the filter outputs are connected.



2) Software Asalab

Before opening the software, check that the amplifier is connected to the computer (USB part of the Fusbi) and that it's switched on.



When the amplifier is switched on, you can see a green LED on it.



a- The USB key for the licence




This USB key (dongle) is necessary to open the software. Don't forget to connect it!

b- How to use Asalab to make an acquisition


Open Asalab.




When you open Asalab you have this window:




Click on Cancel (because you want to make a new acquisition). Click on File and open a New Study.




Now you have this screen:




Now if this is the first time you use this software, you have to create your own Montage. To create your montage click on Montage -> Edit Montage.








This window appears:





Click on New and a new window appears:




Choose a Name (Here I choose ECoG_10E because my montage is for ECoG studies and I have 10 electrodes). Check the box Recording Montage and choose the option Save as Application Montage. In the section Default properties you can choose the Unit, the sensitivity, the colour, the different filters.






When all parameters are chosen. Click on OK.



Now you have this window on your screen:




Now we will see how to create different groups to be able to distinguish the different strips of electrodes. Click on Groups and Click on New.





A new window appears. Choose a name (Strip_4e and Strip_6e for me because I have 2 Strip: one with 4 electrodes and the second with 6 electrodes). You can choose all parameters (colours, filters, sensitivity...)




For each group, we will create electrodes (6 in Strip_6e and 4 in Strip_4e). Select Channels.





Select the first line (= the first electrode) and click on New. Write a label for the electrode and choose the group to apply to this electrode. If you choose default for each parameter (color, sensitivity, filter) the electrode will have all the group parameters that you chose before.




Now we can record signal. To record signal click on Setup Recording.




A new window appears. Choose the parameters. Your new montage appears in the section Montage.









When you chose all the parameters click on OK. To check impedance, Click on impedance Check




When the impedance is good, the colour is green. Now we can record. Choose all the parameters and you click on start :







When your signal is as you can see in the background, you have to change the sensitivity.


Now you have this signal:



To stop Recording click on Stop.




If you want to open your acquisition file, when you will open asalab, in the first window opened, you will select you file:




c- Software parameters

This part is just to explain the main and important options.


Select Tools -> Options...




This window appears:






If you click on each element of the tree (to the left), you will see all the options.









The folder in which the recordings are saved (EEG Folder) must not be modified (or it must be checked that the desired folder has all possible permissions).









3) Issues

a- With Asalab

When we installed the software, the first time, some modules did not install correctly and the problem we faced was that we could not make acquisitions. To solve this problem it was necessary to reinstall the software.

b- With amplifier

We don't know why but when we connect the amplifier by a USB port with the computer, sometimes it is not recognized by the computer. To solve the problem change the USB port.