Scanning protocol (name: INRIA 2018.05.04)

55 slices

resolution 2x2x2mm

freq. dir. L/R

Research → Download

Display → CVs

DELTAx = 0.045 ms

DELTAy = 0.045 ms

DELTAz = 0.045 ms

delta = 0,0322 ms (default setting)

- (1) b-value 400, 6 directions $+ 1b_0$
- (2) b-value 800, 13 directions $+ 1b_0$
- (3) b-value 1550, 29 directions + 1b₀
- (4) b-value 3100, 51 directions $+ 3b_0$

and

(5) T1-weighted

Tips & Tricks

rhimsize = 256

rhrcxres = 256

rhrcyres = 256