

Presurgical MRI data processing

1. Open the following location on NEF:

```
sftp://nef-devel/data/athena/user/monno/connectc/
```

2. Create a directory for the new Patient on NEF. We will call it Patient 1:

```
sftp://nef-devel/data/athena/user/monno/connectc/patient01
```

Note that the directory name can be composed of letters, digits, or the '-' character for a separation. Don't use other characters.

Then, create a subdirectory for storing DICOM files:

```
sftp://nef-devel/data/athena/user/monno/connectc/patient01/dicom
```

3. Copy DICOM images from CD disks to NEF, respecting the structure of directories as in the example below:

```
sftp://nef-devel/data/athena/user/monno/connectc/patient01/dicom/cd1/IMAGES/  
sftp://nef-devel/data/athena/user/monno/connectc/patient01/dicom/cd2/IMAGES/  
sftp://nef-devel/data/athena/user/monno/connectc/patient01/dicom/cd3/IMAGES/
```

Note that the "IMAGES" directory appears in the root directory of every disk recorded by the MRI lab at CHU Nice. It is the only directory that we need.

4. Copy all .sh files from the 'scripts' directory located here:

```
sftp://nef-devel/data/athena/user/monno/connectc/scripts/
```

to our new patient's directory:

```
sftp://nef-devel/data/athena/user/monno/connectc/patient01/
```

5. Login to NEF from the terminal:

```
ssh nef-devel
```

6. Go to our patient's directory:

```
cd /data/athena/user/monno/connectc/patient01
```

7. Run the conversion script:

```
./run_convert.sh
```

Running time is 1-2 days.

8. Use MI-Brain to virtually dissect the interesting fibers. Save them in the .TCK format and place in the following location on NEF:

```
sftp://nef-devel/data/athena/user/monno/connectc/patient01/tracts
```

9. Run the script to create the DICOM images for the neuronavigation system:

```
./run_neuronav.sh
```

Running time is about 15 minutes.

10. Copy the “neuronav” directory to the USB key and load it on the neuronavigation system. The directory can be found on NEF here:

```
sftp://nef-devel/data/athena/user/monno/connectc/patient01/neuronav
```